



CueSee

Hyperbaric CSR

Hyperbaric is transitioning to Hyperbaric CSR

Hyperbaric
Sampling from
Ampoules



Hyperbaric CSR
Sampling with
CueSee® Ringe (CSR)



The Future of Hyperbaric CueSee Has Arrived

**More convenience. More safety.
Same trusted performance.**

For many years, Hyperbaric has provided laboratories with a trusted quality control for verifying blood gas analyzer performance at extremely high pO_2 values. We are excited to announce the next evolution of this trusted product.

Today, Hyperbaric CSR is available alongside the current ampoule version.

Hyperbaric CSR combines the proven Hyperbaric formulation with our innovative CueSee® Ringe sampling technology, making testing easier, safer and more convenient while maintaining the performance laboratories rely on, learn more:*



What changes?

- Same proven Hyperbaric QC formulation
- Integrated CueSee® Ringe sampling system
- Improved ease of use with 15 min open vial stability!
- Safer sample handling
- Reduced risk of handling errors

Transition Timeline

From January 1, 2027

The traditional ampoule version will be discontinued and Hyperbaric CSR will become the standard product available worldwide.

We encourage laboratories to begin evaluating Hyperbaric CSR and transition at a time that best fits their validation and inventory planning.

If you have any questions about the transition or would like assistance implementing Hyperbaric CSR, our team will be happy to help.

Eurotrol

CueSee to the ResCue

*) <https://www.eurotrol.com/en/products/cuesee-hyperbaric-csr/>